


WCS Wheelchair Scale



- **300 kg Capacity**
- **Highly Portable**
- **Simple Operation**
- **Reduces Patient Disturbance**
- **Ergonomically Designed**

AND ...Clearly a Better Value

A&D MERCURY PTY. LTD.

<http://www.andmercury.com.au>





Castor Design for transport

Features

Unique Castor Design

For simple movement of the platform location
(see picture left)

Light Weight Drop Centre Design

Easy access to drop platform surface via entry ramp

Your Choice of Powder Coated or Anodised Finish

Visually appealing for your location

Dual Power Supply

Unit can be operated on AC or battery power

300 kg Capacity - 100 gram Accuracy

Suitable for all wheelchairs

Separate Indicator Stand


Allows maximum flexibility in set up

Keyboard Tare Entry

Subtracts the weight of the wheelchair

Large Weights Area

Suitable for both motorised and non-motorised wheelchairs



Dimensions

WEIGHING AREA HEIGHT 50mm

AND

...Clearly a Better Value

A&D MERCURY PTY. LTD.

ACN 007 556 809

Head Office

32 Dew Street

THEBARTON

South Australia 5031

Telephone (08) 8301 8100

Facsimile (08) 8352 7409

Victorian Office

Unit 4, 313 Arden Street

KENSINGTON

Victoria 3031

Telephone (03) 9372 1522

Facsimile (03) 9372 1193

New South Wales Office

Unit 4, 14 Abbott Road

SEVEN HILLS

New South Wales 2147

Telephone (02) 9674 5466

Facsimile (02) 9674 2544

Email: salesinquiry@andmercury.com.au

Internet page: <http://www.andmercury.com.au>

Your Dealer